
Pharmaceutical Filtration Solutions

Our Membrane For Your Solutions

Membrane Solutions



Enterprise Certification







NSF Certification

CE Certification

FDA Certification







products. It is a high-tech enterprise integrating R&D, production and sales.

Membrane Solutions owns the American patent technology of independent research and development of filter membrane, independently researches and develops and produces Nylon, PVDF, PTFE, PES, etc. (membrane pore size is 0.02-10µm) microporous membrane and membrane filter

INTRODUCTION

Membrane Solutions, founded in 2006, is a


world-renowned supplier of membrane filtration

PTFE, PES, etc. (membrane pore size is 0.02-10µm) microporous membrane and membrane filter cartridge. The corresponding products have obtained ISO9001 and CE certification. And in 2013 passed the US FDA and NSF certification.

The company specializes in the development, production and manufacture of various applied membrane filtration products, such as microfiltration membrane products, membrane pleated filter cartridges, small filters, capsule filters and other products. Membrane Solutions provides professional overall filtration solutions for applications in the global biopharmaceutical field, such as purification and filtration of active pharmaceutical ingredients/reagents, blood products or vaccine products. At the same time, it also provides the most professional fluid filtration and purification solutions for customers in the fields of medical equipment, food and beverage, laboratory analysis, microelectronics industry, etc.



ISO Certification



Registered Trademark Certification

Validation Service Center

As a means of reducing or removing microorganisms (bacteria, mycoplasmas, bacteriophages, endotoxins, etc.) in sterile pharmaceuticals, filtration technology plays a pivotal role in the quality of pharmaceuticals. Membrane Solutions provides products and services that help pharmaceutical manufacturers better comply with industry regulations (GMP) and their stringent standards. The validation support service simplifies the validation process for customers in terms of filter regulations, helping customers reduce system costs.

Membrane Solutions validation center has perfect management procedures, standard validation process and operating specifications, operates in accordance with the company's ISO 9001 quality system requirements, and establishes the management system and technical requirements of the validation center in accordance with CNAS-CL 01:2018 "Accreditation Standard for Testing and Calibration Laboratory Competence". The validation center consists of comprehensive microbial analysis laboratory, comprehensive chemical analysis laboratory, filtration performance laboratory and electron microscope analysis room, which provides complete testing and validation services for food and beverage, pharmaceutical companies and filter manufacturers.









Membrane Solutions Process Validation

Validation Service Scope (Including but not limited to):

- 1.Bacteria Survival Test (Non-final Sterilized Products are required)
- 2.Bacteria Challenge Test (Non-final Sterilized Products are required)
- 3.Product Wettability/Integrity Test
- 4. Chemical Compatibility Test (Non-final Sterilized Products are required)
- 5. Dissolution Test (Non-final Sterilized Products are required)
- 6.Adsorption Test
- 7.Filterability







Membrane Solutions Validation Guide

Chemical Compatibility

Cma-!	Chamical Name	Nylon6,6	PVDF	PTFE	PES	PP	Silicone	EPDM	Viton	DC/D
Species	Chemical Name Glacial Acetic Acid	, .								PS/P\
	Concentrated Hydrochloric Acid	NR NR	R R	R R	NR R	R	LR	R NR	NR NR	R R
	Hydrochloric Acid (6N)	NR	R	R	R	R R	NR NR	NR	R	R
	Concentrated Nitric Acid	NR	R	R	-	R	NR	NR	R	R
Acid	Nitric Acid (6N)	NR	R	R	_	R	LR	NR	R	R
	Concentrated Phosphoric Acid	NR	R	R	-	R	NR	R	R	R
	Concentrated Sulfuric Acid	NR	R	R	NR	R	NR	NR	R	R
	Hydrofluoric Acid (6N)	NR	R	R	-	NR	NR	NR	-	R
	Ammonium Hydroxide (1N)	R	LR	R	R	R	R	R	R	R
	Ammonium Hydroxide (3N)	R	NR	R	R	R	R	LR	LR	R
Alkali	Potassium Hydroxide (3N)	R	LR	R	R	R	LR	R	R	R
	Sodium Hydroxide (3N)	R	LR	R	R	R	R	R	R	R
	Sodium Hydroxide (6N)	R	NR	R	R	R	R	R	R	R
	Amyl Alcohol	R	R	R	R	R	NR	R	R	R
	Benzyl Alcohol	R	R	R	R	R	LR	R	R	R
Alcohols	Butanol	R	R	R	R	R	R	R	LR	R
	Isopropyl Alcohol	R	R	R	R	R	R	R	R	R
	Methanol	LR	R	R	R	R	R	R	NR	R
	Acetone	R	LR	R	NR	R	NR	R	NR	R
Ketone	Cyclohexanone	-	LR	R	NR	R	NR	NR	NR	R
Recone	Methyl Ethyl Ketone	LR	LR	R	-	R	NR	R	NR	R
	Methyl Isobutyl Ketone	LR	LR	R	NR	R	NR	R	NR	R
	Cottonseed Oil	R	R	R	-	R	R	R	R	R
Oil	Lubricating Oil	R	R	R	NR	R	R	R	R	R
	Peanut Oil	R	R	R	-	R	R	R	R	R
	Sesame Oil	R	R	R	R	R	R	R	R	R
Aromatic	Benzene	LR	LR	R	LR	NR	NR	NR	R	R
Hydrocarbon	Toluene	NR	LR	R	NR	NR	NR	NR	R	R
	Xylene Carban Tatranblavida	LR	LR	R	NR	NR	NR	NR	R	R
	Carbon Tetrachloride	LR	LR	R	LR	LR	NR	NR	R	R
	Chloroform Ethylene Dichloride	LR	LR	R	NR	LR	NR	NR	R	R
	Freon TF	LR	LR	R	NR	LR	NR	NR	LR	R
Halogenated Hydrocarbon	Freon TMC	R	R	R R	R	LR	NR	NR	R	R
nyur ocur bon	Methylene Chloride	LR NR	LR LR		NR	LR	NR NR	NR NR	LR LR	R
	Perchloroethylene	NK -	LR	R R	NR LR	LR	NR	NR	R R	R R
	Trichloroethylene	- LR	NR	R	LR	LR LR	NR	NR	R	R
	Ethylene Glycol	R	R	R	LR	R	R	R	R	R
Ethylene	Glycerin	R	R	R	LR	R	R	R	R	R
Glycols	Propylene Glycol	R	R	R	LR	R	R	R	R	R
	Ether	NR	R	R	R	LR	LR	NR	NR	R
February	Isopropyl Ether	-	R	R	-	R	NR	NR	NR	R
Ether	Two 4-dioxane	R	R	R	-	R	NR	NR	NR	R
	Tetrahydrofuran	NR	LR	R	NR	LR	NR	NR	NR	R
	Amyl Acetate	LR	R	R	-	R	NR	R	NR	R
	Butyl Acetate	LR	R	R	-	LR	NR	R	R	R
	Cellulose Acetate	-	R	R	R	R	NR	R	NR	R
Esters	Ethyl Acetate	LR	R	R	LR	LR	R	R	R	R
	Methyl Acetate	LR	R	R	NR	R	-	R	R	R
	Isopropyl Acetate	-	R	R	R	R	LR	R	NR	R
	Aniline	LR	R	R	NR	LR	NR	R	R	R
	Dimethylformamide	R	NR	R	NR	R	R	R	NR	R
	Formaldehyde 37%	R	R	R	R	R	R	R	NR	R
	Gasoline	LR	LR	LR	R	LR	NR	R	R	LR
	Hexane (no water)	-	LR	LR	LR	LR	NR	NR	R	LR
Other	Kerosene	-	R	R	R	R	NR	NR	R	R
Other	Phenol	R	R	R	NR	R	NR	NR	R	R
	Pyridine	LR	R	R	NR	LR	NR	R	NR	R
	Turpentine Oil	-	R	R	R	LR	NR	NR	R	R
	Water	R	R	R	R	R	LR	R	R	R
	Acetonitrile	LR	R	R	R	LR	-	R	NR	R
	Nickel Sulfate Solution	R	R	R		R	R	R	_	R

Note: This table is for reference only



CONTENT




•	Introduction	01
♦	Enterprise Certification	02
•	Validation Service Center	03
♦	Chemical Compatibility Table	04
•	Flow Diagram of Filtration Process in Biopharmaceutical Industry	07
•	PFA All Fluorine Filter Cartridge	09
•	All Fluorine Filter Cartridge with PVDF Cage	10
•	PTFE	11
	FG Series FBG Series FDG Series FL Series	
•	NY	16
	PN Series MZ Series	
•	PP	18
	PPC Series PPD Series PPA Series	
•	ACR Series Activated Carbon Fiber Filter Cartridge	21
•	PVDF	22
	PVL Series PVG Series	

•	PES	24
	SP Series SPB Series SPD Series, SPP Series	
•	Mini Pleated Filter Cartridge	28
•	PP Melt Blown Filter Cartridge	29
•	String Wound Filter Cartridge	30
•	High Flow PP Pleated Filter Cartridge	31
•	Nylon Monofilament Filter Bag	32
•	Needle Felt Filter Bag	33
•	Depth Filter Sheets and Module	34
•	Metal Pleated Filter Cartridge	35
•	Metal Titanium Rod Filter Cartridge	36


Flow Diagram of Filtration Process in Biopharmaceutical Industry

Antibody preparation process Centrifugalization Respirator/ Buffer filtration Culture medium filtration Clarification filtration Pre-column/ pre-ultrafiltration protective filtration Respirator/ gas filtration Buffer filtration Buffer filtration Buffer filtration Pre-column/ pre-ultrafiltration Buffer filtration Concentration & Exchange protective filtration Affinity chromatography Purification Affinity chromatography Virus inactivation Buffer filtration Virus removal pre-filtering Sterilizing filtration Virus filtration Concentration & Exchange


Serological process



Flow Diagram of Filtration Process in Biopharmaceutical Industry



Injection process



PFA All Fluorine Filter Cartridge with Strong Chemical Compatibility

ACG Series

ACG series all fluorine filter cartridge is welded and sealed by ultra-pure PFA resin with high chemical compatibility and high cleanliness PTFE membrane by high-temperature hot-melt technology. The ACG all fluorine filter cartridge does not contain any adhesive. It is washed with $18.2 M\Omega. \text{cm}$ ultra-pure water after production and is available for pre-wetting packaging to ensure very low precipitation levels in harsh use environment, providing good retention efficiency and service life while minimizing the impact on the filtrate.

Technical Parameters


	Filter Membrane	Hydrophobic PTFE			
	Support/Drainage	PFA			
Matariala af	Cage/End Cap	PFA			
Materials of Construction	Core	PFA			
	Adaptor	PFA			
	O-rings	PS/PV			
	Sealed by hot melt welding, no adhesive				
	Outer Diameter	68mm			
Filter	Inner Diameter	33mm			
Dimensions	Length	2.5-40inches			
	Filtration Area	≥0.62m²/10inches			
Operating	Maximum Differential Pressure (Reverse)	4.3 bar @ 25°C,0.5 bar @ 170°C			
Conditions	Maximum Differential Pressure (Forward)	3.0 bar @ 25°C			

Pure Water Flow Rate

- Filtration of strong acid and alkali solutions
- Filtration of strong oxidizing solutions
- Filtration of ultra-pure chemicals
- Filtration of high temperature chemicals
- Sterilization filtration of cleaning solutions and disinfectants



Typical Liquid Flow Rates





Features

- All fluorine structure, strong chemical compatibility
- High cleanliness, pre-wetting packaging is available
- · High flow rate, long service life
- High retention efficiency



Typical Liquid Flow Rates





Features

- Structure with all fluorine materials
- High cleanliness, pre-wetting packaging is available
- High flow rate, long service life
- High retention efficiency

Hydrophobic All Fluorine Filter Cartridge with PVDF Cage

FFV Series

FEV series all fluorine filter cartridge adopts PVDF cage and ECTFE support, and the high porosity PTFE membrane is welded and sealed by high-temperature hot-melt technology. The all fluorine structure provides this product with excellent chemical compatibility and cleanliness. FEVG series can be widely used in the filtration and purification of a variety of strong corrosive solutions. The hydrophilic type can be widely used in the filtration and purification of a wide range of high surface tension and strong corrosive solutions. The hydrophilic filter medium can be used without pre-wetting, which greatly improves the convenience of use.

Technical Parameters

	Filter Membrane	PTFE			
	Support/Drainage	ECTFE			
	Cage/End Cap	PVDF			
Materials of Construction	Core	PVDF			
Construction	Adaptor	PVDF			
	O-rings	PS/PV			
	Sealed by hot melt welding, no adhesive				
	Outer Diameter	68mm			
Filter	Inner Diameter	33mm			
Dimensions	Length	10-40inches			
	Filtration Area	≥0.70m²/10inches			
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 120°C			
Conditions	Maximum Differential Pressure (Forward)	3.4 bar @ 25°C			
Quality Control	ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21				

Typical Applications

- Filtration of strong acid solutions
- Filtration of strong oxidizing solutions
- Filtration of ultra-pure chemicals
- Filtration of high temperature chemicals
- Sterilization filtration of cleaning solutions and disinfectants

• High cleanliness raw materials

Ordering Information

	Length	Pore Size	Adaptor	Seal Material
CRACG	002 = 2.5" 004 = 4" 005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	002 = 0.02μm 003 = 0.03μm 005 = 0.05μm 010 = 0.10μm 020 = 0.20μm 045 = 0.45μm 100 = 1.0μm 500 = 5.0μm H10 = 10.0μm	0=DOE 2=222/Flat 3=222/Fin 6=226/Fin 7=226/Flat	PS=FEP Encapsulated Silicone PV==FEP Encapsulated FPM/FKM

Ordering Information

		Length	Pore Size	Adaptor	Seal Material
CRPTEV	G=Hydrophobic L=Hydrophilic	010 = 10" 020 = 20" 030 = 30" 040 = 40"	005 = 0.05μm 010 = 0.10μm 020 = 0.20μm 045 = 0.45μm 100 = 1.0μm 500 = 5.0μm H10 = 10.0μm	0=DOE 2=222/Flat 3=222/Fin 6=226/Fin 7=226/Flat	PS=FEP Encapsulated Silicone PV==FEP Encapsulated FPM/FKM

• Proven particle retention efficiency and bacteria retention efficiency

Hydrophobic PTFE Pleated Filter Cartridge with Extremely Strong Chemical Compatibility


FG Series

FG series PTFE pleated filter cartridge is made of high porosity hydrophobic PTFE membrane, which has the characteristics of high flux and long service life. Its broad chemical resistance makes the FG series filter cartridge particularly suitable for the filtration of a wide variety of low surface tension chemicals and gases.


Technical Parameters

	Filter Membrane	Hydrophobic PTFE
	Support/Drainage	PP
	Cage/End Cap	PP
Materials of	Core	PP/316L Stainless Steel
Construction	Adaptor	PP/PP+Stainless Steel
	O-rings	S/E/N/F
	Sealed by hot melt wel	ding, no adhesive
	Outer Diameter	69mm
Filter	Inner Diameter	33mm
Dimensions	Length	5-40inches
	Filtration Area	≥0.36m² (5inches), ≥0.74m²/10inches
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C
Biological	Endotoxins	<0.25 EU/mL
Safety	Extractables(WFI)	<15 mg/10 inches
Sterilizable	Steam In-place	140°C 30min differential pressure < 0.3 bar
	Autoclave	121°C 30min
Quality Control	100% integrity testFilter cartridge is fullAll materials complyBacteria (Pseudomo	ty management system ly traceable with relevant requirements of CFR21 nas diminuta) Aerosol Challenge Test Tr≥10 ⁸ for Monodisperse PSL particles≥99.99%

Typical Liquid Flow Rates



Typical Air Flow Rates



Typical Applications

- Filtration of low surface tension solvents
- Filtration of strong acid and alkali liquid
- Filtration of strong oxidizing liquid
- ${}^{\bullet}\:\:$ Sterilization filtration of compressed air and nitrogen
- $\bullet \;\;$ Sterilization filtration of fermentation tank, storage tank, batching tank
- Filtration of solvents and disinfectants

Features `

- $_{\bullet}\,$ Extensive chemical compatibility
- High porosity PTFE membrane
- High flux, low pressure difference
- Very high cost performance


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRPTFEG	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	G10=0.01μm(gas) 010=0.10μm 020=0.20μm 045=0.45μm 100=1.0μm 500=5.0μm H10=10.0μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel


Cost-effective PTFE Pleated Filter Cartridge for GMP Gas Sterilizing Filtration Requirements



Typical Liquid Flow Rates



Typical Air Flow Rates



Typical Applications

- $\bullet \ \ \mbox{Filtration of low surface tension solvents/disinfectants}$
- Sterilization filtration of fermentation tank, storage tank, batching tank
- Filtration of strong acid and alkali/oxidizing liquid
- Sterilization filtration of compressed air/nitrogen

Features \

- · Extensive chemical compatibility
- Proven sterilization efficiency
- High flux, low precipitation
- Water intrusion test can be used

FBG Series

FBG series PTFE pleated filter cartridge is made of naturally hydrophobic sterilizing-grade PTFE membrane, which has excellent uniformity of pore and proven sterilization efficiency, and can provide reliable sterility guarantee for high humidity gas.

Technical Parameters

	Filter Membrane	Hydrophobic PTFE
	Support/Drainage	PP
	Cage/End Cap	PP
Materials of	Core	PP/316L Stainless Steel
Construction	Adaptor	PP/PP+Stainless Steel
	O-rings	S/E/N/F
	Sealed by hot melt w	velding, no adhesive
	Outer Diameter	69mm
Filter	Inner Diameter	33mm
Dimensions	Length	5-40inches
	Filtration Area	≥0.36m² (5inches),≥0.74m²/10inches
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C
Biological	Endotoxins	<0.25 EU/mL
Safety	Extractables	<15 mg/10 inches
Sterilizable	Steam In-place	140°C 30min differential pressure < 0.3 bar
	Autoclave	121°C 30min
Quality Control	100% integrity testFilter cartridge is fuAll materials complBacteria (Pseudome	lity management system Illy traceable y with relevant requirements of CFR21 onas diminuta) Aerosol Challenge Test TR≥10 ⁸ onas diminuta) Liquid Challenge Test TR≥10 ¹¹

	Length	Pore Size	Adaptor	Seal Material	Core
CRPTFEBG	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	G10 = 0.01μm(gas) 010 = 0.10μm 020 = 0.20μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

PTFE Pleated Filter Cartridge for GMP Gas Sterilizing Filtration Requirements


FSG Series

FSG series PTFE pleated filter cartridge is made of strong hydrophobic and uniformly distributed pore size PTFE membrane. Its proven sterilization ability and stable chemical properties enable the filter to provide excellent sterility protection for high humidity gas and low surface tension liquid even after 150 times of moist heat sterilization at 121°C.


Technical Parameters

	Filter Membrane	Hydrophobic PTFE		
	Support/Drainage	PP		
	Cage/End Cap	PP		
Materials of Construction	Core	PP/316L Stainless Steel		
Construction	Adaptor	PP/PP+Stainless Steel		
	O-rings	S/E/N/F		
	Sealed by hot melt welding, no adhesive			
	Outer Diameter	69mm		
Filter	Inner Diameter	33mm		
Dimensions	Length	5-40inches		
	Filtration Area	≥0.36m² (5inches),≥0.74m²/10inches		
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C		
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C		
Biological	Outer Diameter	<0.25 EU/mL		
Safety	Inner Diameter	<15 mg/10 inches		
Sterilizable	Length	140°C 30min differential pressure < 0.3 bar		
Sterilizable	Filtration Area	121°C 30min		
Quality Control	ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Bacteria (Pseudomonas diminuta) Aerosol Challenge Test TR≥10 ⁸ Bacteria Liquid Challenge Test TR≥10 ¹¹			

Typical Liquid Flow Rates



Typical Air Flow Rates



Typical Applications

- Sterilization filtration of compressed air/nitrogen
- Sterilization of fermentation tank, batching tank and breathing filter
- Filtration of aseptic packaging air
- Sterilization filtration of solvents/disinfectants

Features

- Extensive chemical compatibility
- Proven sterilization efficiency
- High reliability
- Water intrusion test can be used


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRPTFESG	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	G10 = 0.01µm(gas) 005 = 0.05µm 010 = 0.10µm 020 = 0.20µm 045 = 0.45µm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=Fluorine Rubber	Blank=PP R=Reinforced PP S=316L Stainless Steel


Double-layer PTFE Pleated Filter Cartridge for GMP Gas Sterilizing Filtration Requirements



Typical Liquid Flow Rates



Typical Air Flow Rates



Typical Applications

- Sterilization filtration of compressed air/nitrogen
- Sterilization of fermentation tank, batching tank and breathing filter
- Filtration of aseptic packaging air
- Sterilization filtration of solvents/disinfectants

Features `

- Extensive chemical compatibility
- High sterilization efficiency
- High flux, low pressure difference
- Water intrusion test can be used

FDG Series

FDG series PTFE pleated filter cartridge is made of double-layer Sterilizing-grade PTFE membrane, which has higher retention efficiency and can meet the requirements of phage retention and other harsh conditions to provide extremely reliable sterility guarantee for gases.

Technical Parameters

	Filter Membrane	Hydrophobic PTFE
	Support/Drainage	PP
	Cage/End Cap	PP
Materials of	Core	PP/316L Stainless Steel
Construction	Adaptor	PP/PP+Stainless Steel
	O-rings	S/E/N/F
	Sealed by hot melt we	elding, no adhesive
	Outer Diameter	69mm
Filter	Inner Diameter	33mm
Dimensions	Length	5-40inches
	Filtration Area	≥0.36m² (5inches),≥0.74m²/10inches
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C
Biological	Endotoxins	<0.25 EU/mL
Safety	Extractables	<15 mg/10 inches
	Steam In-place	140°C 30min differential pressure<0.3 ba
Sterilizable	Autoclave	121°C 30min
Quality Control	100% integrity testFilter cartridge is fullAll materials complyBacteria (Pseudomo	ty management system ly traceable with relevant requirements of CFR21 nas diminuta) Aerosol Challenge Test TR≥10 ¹ nas diminuta) Liquid Challenge Test TR≥10 ¹

	Length	Pore Size	Adaptor	Seal Material	Core
CRPTFEDG	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	G03 = 0.003μm (gas) 010 = 0.10μm 020 = 0.20μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Hydrophilic PTFE Pleated Filter Cartridge with High Flow Rate and Strong Chemical Compatibility

FL Series

FL series PTFE pleated filter cartridge is made of hydrophilic PTFE membrane, which can be used for filtration of various strong acid and alkali and oxidizing solutions without pre-wetting.


Technical Parameters

	Filter Membrane	Hydrophobic PTFE	
	Support/Drainage	PP	
	Cage/End Cap	PP	
Materials of	Core	PP/316L Stainless Steel	
Construction	Adaptor	PP/PP+Stainless Steel	
	O-rings	S/E/N/F	
	Sealed by hot melt we	elding, no adhesive	
	Outer Diameter	69mm	
Filter	Inner Diameter	33mm	
Dimensions	Length	5-40inches	
	Filtration Area	≥0.36m² (5inches),≥0.74m²/10inches	
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C	
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C	
Biological	Endotoxins	<0.25 EU/mL	
Safety	Extractables	<15 mg/10 inches	
6: 11: 11	Steam In-place	140°C 30min differential pressure < 0.3 bar	
Sterilizable	Autoclave	121°C 30min	
Quality Control	 ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Retention Efficiency for Monodisperse PSL particles≥99.99% 		

15



Typical Liquid Flow Rates



Features

- Extensive chemical compatibility
- Good corrosion resistance, oxidation resistance, heat resistance
- High porosity PTFE membrane, no pre-wetting required
- High flux, low precipitation, low pressure difference
- 100% integrity test ensures the effect of sterilization

Typical Applications

- Filtration of high temperature and high viscosity solvents
- Filtration of strong acid liquid
- Filtration of strong oxidizing liquid
- Filtration of high-purity chemicals
- Filtration of LVP
- Filtration of biological reagents

Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRPTFEL	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	010 = 0.10μm 020 = 0.20μm 045 = 0.45μm 100 = 1.0μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Nylon66 Pleated Filter Cartridge with High Strength and High Temperature Resistance




PN Series

PN Series Nylon pleated filter cartridge is made of hydrophilic nylon66 membrane with lining, which has excellent chemical compatibility and mechanical strength. After special cleaning steps, the filter membrane has low precipitation and no shedding. It can be widely used in the filtration of strong alkali, weak acid and organic solvents.

Technical Parameters

	Filter Membrane	Nylon66PP
	Support/Drainage	PP
	Cage/End Cap	PP/316L Stainless Steel
Materials of	Core	PP/PP+Stainless Steel
Construction	Adaptor	S/E/N/F
	O-rings	
	Sealed by hot melt we	elding, no adhesive
	Outer Diameter	69mm
Filter	Inner Diameter	33mm
Dimensions	Length	5-40inches
	Filtration Area	≥0.33m² (5inches),≥0.62m²/10inches
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C
Biological	Endotoxins	<0.25 EU/mL
Safety	Extractables	<40 mg/10 inches
	Steam In-place	121°C 30min differential pressure<0.3 bar
Sterilizable	Autoclave	121°C 30min
Quality Control	100% integrity testFilter cartridge is fulAll materials complyBacteria (Pseudomo	ity management system lly traceable y with relevant requirements of CFR21 onas diminuta) Liquid Challenge Test TR≥10 ¹¹ for Monodisperse PSL particles≥99.99%

Typical Liquid Flow Rates



Features

- Naturally hydrophilic, easy to wet
- Extensive chemical compatibility, resistant to organic solvents • Good integrity to ensure the filter sterile-grade particle retention
- Low pressure difference, high flux, long service life

Typical Applications

- Sterilization filtration of sterile APIs, antibiotic solutions, etc.
- Filtration of alcohol, soft drinks, pure water, etc.
- Filtration of high- purity chemicals and organic solvents

• Sterilization filtration of physiological saline and reagent water • Sterilization filtration of SVP and LVP

	Length	Pore Size	Adaptor	Seal Material	Core
CRNY	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	02X = 0.20μm + 0.20μm 020 = 0.20μm 045 = 0.45μm 065 = 0.65μm 100 = 1.0μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS inser	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Positively Charged Modified Nylon 66 Pleated Filter Cartridge

MZ Series

The positively charged modified nylon membrane used in the MZ series has a stable positive charge point in water-based solutions. The positive charge of the nylon membrane loaded with high specific surface area can efficiently adsorb negatively charged impurities, whose size is much smaller than the pore size of the filter


Technical Parameters

	Filter Membrane	Nylon66PP
	Support/Drainage	PP
	Cage/End Cap	PP/316L Stainless Steel
Materials of	Core	PP/PP+Stainless Steel
Construction	Adaptor	S/E/N/F
	O-rings	
	Sealed by hot melt we	elding, no adhesive
	Outer Diameter	69mm
Filter	Inner Diameter	33mm
Dimensions	Length	5-40inches
	Filtration Area	≥0.33m² (5inches),≥0.62m²/10inches
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C
Biological	Endotoxins	<0.25 EU/mL
Safety	Extractables	<40 mg/10 inches
	Steam In-place	121°C 30min differential pressure < 0.3 bar
Sterilizable	Autoclave	121°C 30min
Quality Control	100% integrity testFilter cartridge is fulAll materials complyBacteria Liquid Cha	with relevant requirements of CFR21

17



Typical Liquid Flow Rates



Features

- Highly hydrophilic, easy to wet
- Positive charge on membrane surface
- Good integrity to ensure the filter sterile-grade particle retention efficiency
- No fiber shedding

Typical Applications

- Sterilization filtration of sterile APIs, antibiotic solutions, etc.
- Sterilization filtration of intermediates, soft drinks, pure water, etc.
- Filtration of high- purity chemicals and organic solvents
- Sterilization filtration of physiological saline, reagents and water for injection
- Sterilization filtration of SVP and LVP


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRNYZ	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	005 = 0.05μm· 010 = 0.10μm 020 = 0.20μm 045 = 0.45μm 100 = 1.0μm 500 = 5.0μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=Fluorine Rubber	Blank=PP R=Reinforced PP S=316L Stainless Steel

PP Pleated Filter Cartridge with High Flux and High Flow Rate



Typical Liquid Flow Rates



Features `

- Extensive chemical compatibility
- High efficiency, high flux, low pressure difference, long service life
- High efficiency particle retention rate to protect terminal filter sterilization safety
- No fiber shedding

Typical Applications

- Filtration of colloidal products
- Filtration of fermentation broth and high viscosity materials
- Prefiltration of culture media and blood products
- Filtration of chemical intermediates, organic solvents, pure water

PPC Series

PPC series PP pleated filter cartridge is made of high fluffy polypropylene filter material, which has the characteristics of high flow rate and high dirt holding capacity. It has excellent cost performance as clarification filter and pre-filter. The full PP structure enables the filter cartridge to be applied to the filtration of various acid and alkali and organic solvents.

Technical Parameters

Filter Membrane Support/Drainage PP Cage/End Cap PPP Core PP/316L Stainless Steel Adaptor PP/PP+Stainless Steel O-rings S/E/N/F Sealed by hot melt welding, no adhesive Outer Diameter Dimensions Inner Diameter S-40inches Filtration Area Perssure (Reverse) Maximum Differential Pressure (Reverse) Maximum Differential Safety Endotoxins Sterilizable Sterilizable Autoclave Filter Conditions Filtration Area Support/Drainage PP				
Cage/End Cap PP Core PP/316L Stainless Steel Adaptor PP/PP+Stainless Steel O-rings S/E/N/F Sealed by hot melt welding, no adhesive Outer Diameter 69mm Inner Diameter 33mm Length 5-40inches Filtration Area ≥0.25m² (5inches), ≥0.52m²/10inches Filtration Area ≥0.25m² (5inches), ≥0.52m²/10inches Operating Conditions Maximum Differential Pressure (Reverse) Maximum Differential Pressure (Forward) Biological Safety Extractables <0.25 EU/mL Sterilizable Autoclave 121°C 30min differential pressure <0.3 bar and the pressure <0.3		Filter Membrane	PP	
Materials of Core PP/316L Stainless Steel Construction Core PP/PP+Stainless Steel O-rings S/E/N/F Sealed by hot melt welding, no adhesive Outer Diameter 69mm Inner Diameter 33mm Length 5-40inches Filtration Area ≥0.25m² (5inches),≥0.52m²/10inches Operating Conditions Maximum Differential Pressure (Reverse) Maximum Differential Pressure (Forward) 2.1 bar @ 25°C Biological Safety Endotoxins <0.25 EU/mL		Support/Drainage	PP	
Core PP/316L Stainless Steel Adaptor PP/PP+Stainless Steel O-rings S/E/N/F Sealed by hot melt welding, no adhesive Outer Diameter 69mm Inner Diameter 33mm Length 5-40inches Filtration Area ≥0.25m² (5inches),≥0.52m²/10inches Filtration Area ≥0.25m² (5inches),≥0.52m²/10inches Operating Pressure (Reverse) Conditions Maximum Differential Pressure (Reverse) Maximum Differential Pressure (Forward) Safety Endotoxins <0.25 EU/mL Safety Extractables <40 mg/10 inches Sterilizable Autoclave 121°C 30min PP/PP+Stainless Steel PP/PP+Stainless Steel S-/E/N/F Sealed by hot melt welding, no adhesive 5.2 bar © 2.1 bar © 25°C,1.9 bar © 82°C 2.1 bar © 25°C 2.1 bar © 25°C 121°C 30min differential pressure <0.3 bar on the place of the		Cage/End Cap	PP	
Adaptor PP/PP+Stainless Steel O-rings S/E/N/F Sealed by hot melt welding, no adhesive Outer Diameter 69mm Inner Diameter 33mm Length 5-40inches Filtration Area ≥0.25m² (5inches),≥0.52m²/10inches Operating Pressure (Reverse) Conditions Maximum Differential Pressure (Reverse) Conditions Maximum Differential Pressure (Forward) Biological Safety Extractables <0.25 EU/mL Sterilizable Autoclave 121°C 30min Pressure (Soundity management system Filter cartridge is fully traceable	materials of	Core	PP/316L Stainless Steel	
Sealed by hot melt welding, no adhesive Outer Diameter 69mm Inner Diameter 33mm Length 5-40inches Filtration Area ≥0.25m² (5inches),≥0.52m²/10inches Operating Conditions Maximum Differential Pressure (Reverse) 5.2 bar @ 25°C,1.9 bar @ 82°C Maximum Differential Pressure (Forward) 2.1 bar @ 25°C Biological Endotoxins <0.25 EU/mL Safety Extractables <40 mg/10 inches Sterilizable Autoclave 121°C 30min Ouality • ISO 9001:2015 quality management system • Filter cartridge is fully traceable	Construction	Adaptor	PP/PP+Stainless Steel	
Outer Diameter 69mm Inner Diameter 33mm Length 5-40inches Filtration Area ≥0.25m² (5inches),≥0.52m²/10inches Operating Conditions Maximum Differential Pressure (Reverse) 2.1 bar @ 25°C Biological Pressure (Forward) 2.1 bar @ 25°C Endotoxins <0.25 EU/mL Safety Extractables <40 mg/10 inches Sterilizable Autoclave 121°C 30min Outer Diameter 69mm 33mm 5-40inches 5.2 bar @ 25°C,1.9 bar @ 82°C 2.1 bar @ 25°C 40 mg/10 inches 121°C 30min differential pressure <0.3 bar and the pressure <0.3 bar and t		O-rings	S/E/N/F	
Filter Dimensions Length S-40inches Filtration Area Dimensions Deparating Conditions Maximum Differential Pressure (Reverse) Maximum Differential Pressure (Forward) Maximum Differential Pressure (Forward) Maximum Differential Pressure (Forward) Endotoxins Safety Extractables Steam In-place Autoclave Pilter cartridge is fully traceable Samm 5-40inches 5.2 bar @ 25°C, 1.9 bar @ 82°C 2.1 bar @ 25°C 2.1 bar @ 25°C 2.1 bar @ 25°C 2.2 bar @ 25°C 2.3 bar @ 25°C 2.1 bar @ 25°C 2.1 bar @ 25°C 2.2 bar @ 25°C 2.3 bar @ 25°C 2.4 bar @ 25°C 2.5 EU/mL 2.7 c 30min differential pressure < 0.3 bar @ 25°C 2.8 bar @ 25°C 2.9 bar @ 25°C 2.1 bar @ 25°C 2.1 bar @ 25°C 2.1 bar @ 25°C 2.2 bar @ 25°C 2.3 bar @ 25°C 2.4 bar @ 25°C 2.5 EU/mL 3.5 EV/mL 4.6 bar @ 25°C 4.7 bar @ 25°C 2.8 bar @ 25°C 2.9 bar @ 25°C 2.1 b		Sealed by hot melt we	elding, no adhesive	
Dimensions Length 5-40inches Filtration Area Dimensions Length Filtration Area Solutions Maximum Differential Pressure (Reverse) Maximum Differential Pressure (Forward) Biological Safety Endotoxins Sterilizable Sterilizable Autoclave Pilter cartridge is fully traceable S-40inches 5-40inches 5-40inches 5-2 bar @ 25°C,1.9 bar @ 82°C 2.1 bar @ 25°C 2.1 bar @ 25°C 2.1 bar @ 25°C 2.25 EU/mL 40 mg/10 inches 121°C 30min		Outer Diameter	69mm	
Filtration Area ≥0.25m² (5inches),≥0.52m²/10inches Operating Conditions Maximum Differential Pressure (Reverse) Maximum Differential Pressure (Forward) Biological Endotoxins <0.25 EU/mL Extractables <40 mg/10 inches Sterilizable Autoclave 121°C 30min • ISO 9001:2015 quality management system • Filter cartridge is fully traceable	Filter	Inner Diameter	33mm	
Operating Pressure (Reverse) Conditions Maximum Differential Pressure (Reverse) Maximum Differential Pressure (Forward) Biological Endotoxins Safety Extractables Sterilizable Autoclave Ouality Autoclave Sterilizable Maximum Differential Pressure (2.1 bar @ 25°C) 2.1 bar @ 25°C 2.2 bar @ 25°C 2.1 bar @ 25°C 2.2 bar @ 25°C 2.3 bar @ 25°C 2.4 bar @ 25°C 2.5 EU/mL Sterilizable 40 mg/10 inches 121°C 30min 121°C 30min ISO 9001:2015 quality management system Filter cartridge is fully traceable	Dimensions	Length	5-40inches	
Operating Conditions Pressure (Reverse) Maximum Differential Pressure (Forward) Biological Safety Endotoxins Sterilizable Sterilizable Pressure (Reverse) 2.1 bar @ 25°C 2.1 bar @ 25°C 2.25 EU/mL 40 mg/10 inches 121°C 30min differential pressure < 0.3 bar and 25°C Autoclave 121°C 30min ISO 9001:2015 quality management system Filter cartridge is fully traceable		Filtration Area	≥0.25m² (5inches),≥0.52m²/10inches	
Pressure (Forward) 2.1 bar @ 25°C	Operating		5.2 bar @ 25°C,1.9 bar @ 82°C	
Safety Extractables Steam In-place Autoclave 121°C 30min 121°C 30min 121°C 30min • ISO 9001:2015 quality management system • Filter cartridge is fully traceable	Conditions		2.1 bar @ 25°C	
Steam In-place 121°C 30min differential pressure < 0.3 bath and the state of the st	Biological	Endotoxins	<0.25 EU/mL	
Sterilizable Autoclave 121°C 30min ISO 9001:2015 quality management system Filter cartridge is fully traceable	Safety	Extractables	<40 mg/10 inches	
Autoclave 121°C 30min • ISO 9001:2015 quality management system • Filter cartridge is fully traceable	6. 11. 11	Steam In-place	121°C 30min differential pressure<0.3 ba	
Filter cartridge is fully traceable	Sterilizable	Autoclave	121°C 30min	
		Filter cartridge is ful	lly traceable	

	Length	Pore Size	Adaptor	Seal Material	Core
CRPP	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	010 = 0.10µm 020 = 0.20µm 045 = 0.45µm 100 = 1.0µm 300 = 3.0µm 500 = 5.0µm H10 = 10µm H20 = 20µm H30 = 30µm H50 = 50µm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Multi-layer PP Pleated Filter Cartridge with Long Service Life

PPD Series


PPD adopts multi-layer gradient PP fiber material, which combines the advantages of depth filtration and pleated filtration, so that the filter cartridge has high dirt holding capacity and high retention efficiency. The full PP structure makes PPD series pleated filter cartridge suitable for clarification and filtration of various high-load chemicals.

Technical Parameters

	Filter Membrane	PP
	Support/Drainage	PP
	Cage/End Cap	PP
Materials of	Core	PP/316L Stainless Steel
Construction	Adaptor	PP/PP+Stainless Steel
	O-rings	S/E/N/F
	Sealed by hot melt we	elding, no adhesive
	Outer Diameter	69mm
Filter	Inner Diameter	33mm
Dimensions	Length	5-40inches
	Filtration Area	≥0.15m² (5inches),≥0.30m²/10inches
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C
Biological	Endotoxins	<0.25 EU/mL
Safety	Extractables	<40 mg/10 inches
	Steam In-place	121°C 30min differential pressure < 0.3 bar
Sterilizable	Autoclave	121°C 30min
Quality Control	ISO 9001:2015 quality management system Filter cartridge is fully traceable All materials comply with relevant requirements of CF	



Typical Liquid Flow Rates



Features

- Multi-layer pore size gradient filter medium, extensive chemical compatibility
- Combination of depth filter and pleated filter
- High efficiency, high flux, low pressure difference, long service life
- No fiber shedding

Typical Applications

- Filtration of glue, dyestuff, grinding fluid, high-load filtrate
- Filtration of colloidal products
- Filtration of fermentation broth and high viscosity materials
- Prefiltration of culture media, blood products


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRPPD	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	050 = 0.50μm 080 = 0.80μm 100 = 1.0μm 150 = 1.5μm 300 = 3.0μm 500 = 5.0μm H10 = 10μm H15 = 15μm H20 = 20μm H30 = 30μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

PP Pleated Filter Cartridge with Higher Removal Efficiency



Typical Liquid Flow Rates



Features

- Extensive chemical compatibility
- Gradient pore size design
- High efficiency particle retention rate to protect terminal filter sterilization safety
- No fiber shedding

Typical Applications

- Filtration of colloidal products
- Filtration of fermentation broth and high viscosity materials
- Prefiltration of culture media, blood products
- Filtration of chemical intermediates, organic solvents, pure water

PPA Series

PPA series PP pleated filter cartridge is made of multi-layer PP microfiber. With gradient pore size and high strength support, the filter cartridge has high retention efficiency and high flow rate. The full PP structure enables the filter cartridge to be applied to the filtration of various acid and alkali and organic solvents.

Technical Parameters

	Filter Membrane	PP	
	Support/Drainage	PP	
	Cage/End Cap	PP	
Materials of	Core	PP/316L Stainless Steel	
Construction	Adaptor	PP/PP+Stainless Steel	
	O-rings	S/E/N/F	
	Sealed by hot melt we	elding, no adhesive	
	Outer Diameter	69mm	
Filter	Inner Diameter	33mm	
Dimensions	Length	5-40inches	
	Filtration Area	≥0.25m² (5inches),≥0.52m²/10inches	
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C	
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C	
Biological	Endotoxins	<0.25 EU/mL	
Safety	Extractables	<40 mg/10 inches	
	Steam In-place	121°C 30min differential pressure < 0.3 bar	
Sterilizable	Autoclave	121°C 30min	
Quality Control	ISO 9001:2015 quality management system Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Proven microbial and particle retention efficiency		

	Length	Pore Size	Adaptor	Seal Material	Core
CRPPA	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	020 = 0.20µm 045 = 0.45µm 065 = 0.65µm 100 = 1.0µm 300 = 3.0µm 500 = 5.0µm H10 = 10µm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Avtivated Carbon Fiber Filter Cartridge

ACR Series

The food-grade activated carbon fiber used in the ACR series activated carbon fiber winding filter cartridge has extremely high surface area and high porosity. It has the characteristics of fast adsorption rate and high adsorption capacity and can provide ideal retention efficiency.

Technical Parameters

	Filter Membrane	Activated Carbon Fiber		
	Support/Drainage	PP		
	Cage/End Cap	PP		
Materials of	Core	PP/316L Stainless Steel		
Construction	Adaptor	PP/PP+Stainless Steel		
	O-rings	S/E/N/F		
	Sealed by hot melt welding, no adhesive			
	Outer Diameter	69 mm		
Filter Dimensions	Inner Diameter	33 mm		
	Length	5-40inches		
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25℃,1.9 bar @ 82℃		
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C		
Quality Control	ISO 9001:2015 quality management system Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21			

21



Features

- · Activated carbon with high iodine value
- Fast adsorption
- Depth filter
- · Long service life

Typical Applications

- Decolorization of material liquid
- Deodorization of material liquid
- · Clarification and filtration
- Dechlorination of material liquid


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRACR	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	050 = 0.50μm 100 = 1.0μm 150 = 1.5μm 300 = 3.0μm 500 = 5.0μm H10 = 10μm H15 = 15μm H30 = 30μm H50 = 50μm T10 = 100μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Hydrophilic PVDF Pleated Filter Cartridge with Broad Chemical Compatibility

PVL Series

PVL series PVDF pleated filter cartridge is made of unlined hydrophilic PVDF membrane. The filter membrane has the characteristics of symmetrical pore type, uniform hydrophilicity, high porosity and extremely low protein binding, which make the filter cartridge easy to use without pre-wetting. The filter cartridge is suitable for the sterilization filtration of various biological products.



Features

- Low precipitation and have extensive chemical compatibility
- High porosity, very low protein binding
- 100% integrity test ensures the effect of germ removal
- No pre-wetting required

Typical Applications

- Filtration of high-purity chemicals, chemical intermediates,
- strong acid and oxidizing liquid Filtration of SVP and LVP
- Filtration of sterile APIs and ophthalmic preparations
- Filtration of vaccines, biological products, blood products and other high protein content materials

Technical Parameters

	Filter Membrane	Hydrophilic PVDF	
	Support/Drainage	PP	
	Cage/End Cap	PP/316L Stainless Steel	
Materials of	Core	PP/PP+Stainless Steel	
Construction	Adaptor	S/E/N/F	
	O-rings		
	Sealed by hot melt we	elding, no adhesive	
	Outer Diameter	69mm	
Filter	Inner Diameter	33mm	
Dimensions	Length	5-40inches	
	Filtration Area	≥0.32m² (5inches),≥0.65m²/10inches	
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C	
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C	
Biological	Endotoxins	<0.25 EU/mL	
Safety	Extractables	<40 mg/10 inches	
6. 11. 1.1	Steam In-place	121°C 30min differential pressure < 0.5 bar	
Sterilizable	Autoclave	121°C 30min	
Quality Control	 ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Proven bacteria retention efficiency 		

	Length	Pore Size	Adaptor	Seal Material	Core
CRPVDFL	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	010 = 0.10µm 020 = 0.20µm 045 = 0.45µm 065 = 0.65µm 100 = 1.0µm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Hydrophobic PVDF Pleated Filter Cartridge with High Flow Rate and Long Service Life

PVG Series


PVG series PVDF pleated filter cartridge is made of unlined hydrophobic PVDF membrane. The filter membrane has the characteristics of symmetrical pore type, high porosity and extremely low protein binding, which make the filter cartridge has high initial flow rate and longer service life. The filter cartridge is suitable for the filtration of various chemical solvents and biological products.

Technical Parameters

	Filter Membrane	Hydrophobic PVDF	
	Support/Drainage	PP	
	Cage/End Cap	PP	
Materials of	Core	PP/316L Stainless Steel	
Construction	Adaptor	PP/PP+Stainless Steel	
	O-rings	S/E/N/F	
	Sealed by hot melt we	elding, no adhesive	
	Outer Diameter	69mm	
Filter	Inner Diameter	33mm	
Dimensions	Length	5-40inches	
	Filtration Area	≥0.32m² (5inches),≥0.65m²/10inches	
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C	
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C	
Biological	Endotoxins	<0.25 EU/mL	
Safety	Extractables	<40 mg/10 inches	
	Steam In-place	121°C 30min differential pressure<0.5 bar	
Sterilizable	Autoclave	121°C 30min	
Quality Control	 ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Proven bacteria retention efficiency 		



Typical Liquid Flow Rates



Features

- Low precipitation, very low protein binding
- Extensive chemical compatibility
- 100% integrity test ensures the effect of germ removal
- High porosity, long service life

Typical Applications

- Filtration of high-purity chemicals, chemical intermediates, strong acid and oxidizing liquid
- Filtration of SVP and LVP
- Filtration of sterile APIs and ophthalmic preparations
- Filtration of vaccines, biological products, blood products and other high protein content materials


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRPVDFG	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	010 = 0.10μm 020 = 0.20μm 045 = 0.45μm 065 = 0.65μm 100 = 1.0μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

PES Pleated Filter Cartridge with High Flow Rate and Long Service Life



Typical Liquid Flow Rates



Features \

- · Strong chemical compatibility
- High flux, low protein binding
- 100% integrity test
- Long service life

Typical Applications

- Filtration of water for injection, cleaning solutions, purified water, deionized water
- Filtration of LVP, APIs, buffers, eye drops, disinfectants, etc.
- Filtration of vaccines, serums, biological products, antibiotic water-based liquid, etc.
- $\bullet \;\;$ Filtration of red wine, beer, juice, purified water, etc.

SP Series

SP series PES pleated filter cartridge is made of naturally hydrophilic highly asymmetric PES membrane. It has good heat resistance and can withstand hot water sterilization and steam sterilization for a long time. The removal rate of the filter cartridge to typical bacteria in industry passes bath validation. Longer flux can be provided while effectively reducing the load.

Technical Parameters

	Filter Membrane	PES	
	Support/Drainage	PP	
	Cage/End Cap	PP	
Materials of	Core	PP/316L Stainless Steel	
Construction	Adaptor	PP/PP+Stainless Steel	
	O-rings	S/E/N/F	
	Sealed by hot melt we	elding, no adhesive	
	Outer Diameter	69mm	
Filter	Inner Diameter	33mm	
Dimensions	Length	5-40inches	
	Filtration Area	≥0.32m² (5inches),≥0.65m²/10inches	
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C	
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C	
Biological	Endotoxins	<0.25 EU/mL	
Safety	Extractables	<40 mg/10 inches	
Constitution in the	Steam In-place	121°C 30min differential pressure < 0.5 bar	
Sterilizable	Autoclave	121°C 30min	
Quality Control	 ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Retention Efficiency for Monodisperse PSL particles≥99.99% 		

	Length	Pore Size	Adaptor	Seal Material	Core
CRPES	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	002 = 0.02μm 005 = 0.05μm 010 = 0.10μm 020 = 0.20μm 045 = 0.45μm 065 = 0.65μm 080 = 0.80μm 120 = 1.20μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Sterilizing-grade PES Pleated Filter Cartridge with High Flow Rate and High Efficiency

SPB Series


SPB series PES pleated filter cartridge is made of imported asymmetric PES membrane, which has ultra-low protein binding and dissolution/precipitation level. Its good heat resistance enables the filter cartridge to withstand multiple moist heat sterilization. Batch bacteria challenge guarantees the stable bacteria removal ability of the filter cartridge.

Technical Parameters

	Filter Membrane	PES	
	Support/Drainage	PP	
	Cage/End Cap	PP/316L Stainless Steel	
Materials of	Core	PP/PP+Stainless Steel	
Construction	Adaptor	S/E/N/F	
	O-rings		
	Sealed by hot melt we	elding, no adhesive	
	Outer Diameter	69mm	
Filter	Inner Diameter	33mm	
Dimensions	Length	5-40inches	
	Filtration Area	≥0.32m² (5inches),≥0.65m²/10inches	
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C	
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C	
Biological	Endotoxins	<0.25 EU/mL	
Safety	Extractables	<40 mg/10 inches	
	Steam In-place	121°C 30min differential pressure < 0.5 bar	
Sterilizable	Autoclave	121°C 30min	
Quality Control	 ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Bacteria Liquid Challenge Test T_R≥10¹¹ 		



Typical Liquid Flow Rates



Features

- · Strong chemical compatibility
- · High flux, low protein binding
- Proven sterile-grade filter cartridge
- Ultra-low dissolved/precipitated materials

Typical Applications

- Sterilization filtration of water for injection, cleaning solutions, purified water, deionized water, etc.
- Sterilization filtration of LVP, APIs, buffers, eye drops, disinfectants, etc.
- Sterilization filtration of vaccines, serums, biological products, antibiotic water-based liquid, etc.
- Filtration of red wine, beer, juice, purified water, etc.


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRPESB	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	010 = 0.10μm 020 = 0.20μm 045 = 0.45μm 065 = 0.65μm 080 = 0.80μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel

Sterilizing-grade Double-layer PES Pleated Filter Cartridge with Higher Removal Efficiency



Typical Liquid Flow Rates



Features

- Strong chemical compatibility
- · High flux, low protein binding
- Proven sterile-grade filter cartridge
- Higher protection ability of sterilization

Typical Applications

- Sterilization filtration of water for injection, cleaning solutions, purified water, deionized water
- · Sterilization filtration of LVP, APIs, buffers, eye drops, disinfectants, etc.
- Sterilization filtration of vaccines, serums, biological products, antibiotic water-based liquid, etc.
- Filtration of red wine, beer, juice, purified water

SPD series PES pleated filter cartridge is made of double-layer imported asymmetric PES membrane. The double-layer Sterilizing-grade membrane provides a higher level of sterilization protection. It is suitable for harsh sterile terminal filtration.

Technical Parameters

	Filter Membrane	PES				
	Support/Drainage	PP				
	Cage/End Cap	PP/316L Stainless Steel				
Materials of	Core	PP/PP+Stainless Steel				
Construction	Adaptor	S/E/N/F				
	O-rings					
	Sealed by hot melt welding, no adhesive					
	Outer Diameter	69mm				
Filter	Inner Diameter	33mm				
Dimensions	Length	5-40inches				
	Filtration Area	≥0.26m² (5inches),≥0.55m²/10inches				
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C				
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C				
Biological	Endotoxins	<0.25 EU/mL				
Safety	Extractables	<40 mg/10 inches				
	Steam In-place	121°C 30min differential pressure<0.5 ba				
Sterilizable	Autoclave 121°C 30min					
Quality Control	 ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Bacteria Liquid Challenge Test TR≥10¹¹ 					

	Length	Pore Size	Adaptor	Seal Material	Core
CRPESD	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	010 = 0.10μm + 0.10μm 020 = 0.20μm+ 0.20μm 045 = 0.45μm+ 0.45μm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Stee

Sterilizing-grade Double-layer PES Pleated Filter Cartridge with Long Service Life and High Efficiency

SPP Series


SPP series PES pleated filter cartridge is made of double-layer imported asymmetric PES membrane. The double-layer filter membrane is specially designed so that the filter cartridge has reliable sterilization ability and can provide longer flux.

Technical Parameters

	Filter Membrane	PES		
	Support/Drainage	PP		
	Cage/End Cap	PP/316L Stainless Steel		
Materials of	Core	PP/PP+Stainless Steel		
Construction	Adaptor	S/E/N/F		
	O-rings			
	Sealed by hot melt we	elding, no adhesive		
	Outer Diameter	69mm		
Filter	Inner Diameter	33mm		
Dimensions	Length	5-40inches		
	Filtration Area	≥0.26m² (5inches),≥0.55m²/10inches		
Operating	Maximum Differential Pressure (Reverse)	5.2 bar @ 25°C,1.9 bar @ 82°C		
Conditions	Maximum Differential Pressure (Forward)	2.1 bar @ 25°C		
Biological	Endotoxins	<0.25 EU/mL		
Safety	Extractables	<40 mg/10 inches		
6. 11. 11	Steam In-place	121°C 30min differential pressure < 0.3 ba		
Sterilizable	Autoclave	121°C 30min		
Quality Control	 ISO 9001:2015 quality management system 100% integrity test Filter cartridge is fully traceable All materials comply with relevant requirements of CFR21 Bacteria Liquid Challenge Test T_R≥10¹¹ 			



Typical Liquid Flow Rates



Features

- Strong chemical compatibility
- · High flux, low protein binding
- Proven sterile-grade filter cartridge
- Super long service life

Typical Applications

- Sterilization filtration of water for injection, cleaning solutions, purified water, deionized water
- Sterilization filtration of LVP, APIs, buffers, eye drops, disinfectants, etc.
- Sterilization filtration of vaccines, serums, biological products, antibiotic water-based liquid, etc.
- Filtration of red wine, beer, juice, purified water


Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Core
CRPESP	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	010=0.10µm+0.20µm 020=0.20µm+0.80µm 045=0.45µm+1.0µm	0=DOE 1=222/Flat with 316L SS insert 2=222/Flat 3=222/Fin 4=222/Fin with 316L SS insert 5=226/Fin with 316L SS insert 6=226/Fin 7=226/Flat 8=226/Flat with 316L SS insert	S=Silicone E=EPDM N=Nitrile F=FKM	Blank=PP R=Reinforced PP S=316L Stainless Steel


Mini Pleated Filter Cartridge


JPF Series

JPF series mini pleated filter cartridges are used for laboratory research, small-scale process validation in the pharmaceutical industry and other sterilization filtration, particle control, pre-filtration, etc. We can provide microporous membrane with different materials and removal ratings and product selection of different installation interface to meet various needs of customers.



Typical Liquid Flow Rates





Technical Parameters

	Filter Membrane	Hydrophobic PTFE/PP/PES		
	Support/Drainage	PP		
	Cage/End Cap	PP		
Materials of	Core	PP		
Construction	Adaptor	PP		
	O-rings	S/E/F		
	Sealed by hot melt welding, no adhesive			
	Outer Diameter	J25/J50: 56mm F25:42mm		
Filter	Length	2.5inches, 5inches		
Dimensions	Filtration Area	J25≥0.1m²,J50≥0.2m²,F25≥0.045m²		
	Maximum Differential Pressure (Forward)	4.2 bar @ 25°C,1.5 bar @ 82°C, 0.5 bar @ 95°C		
Operating Conditions	Maximum Differential Pressure (Reverse)	2.1 bar @ 25°C		
	Steam In-place	121°C 30min differential pressure < 0.3 bar		
	Autoclave	121°C 30min		

Features

- Interpolation is compact and easy to use
- Proven sterilization efficiency
- High flow rate, long service life

Typical Applications

- · Sterilization filtration of compressed air, nitrogen, etc.
- Filtration of aseptic packaging air
- Withstand multiple moist heat sterilization Sterilization filtration of solvents, disinfectants, etc.
 - Sterilization filtration of biological preparations, culture medium, sterile APIs, antibiotic liquid, etc.

	Filter Membrane	Length	Pore Size	Adaptor	Seal Material
CR	PP PTFEG=Hydrophobic PTFE PTFESG=Sterile-grade Hydrophobic PTFE PES PESS PESB=Sterile-grade PES	F25=2.5" J25=2.5" J50=5"	010=0.10µm 020=0.20µm 045=0.45µm 100=1.00µm 300=3.00µm 500=5.00µm H10=10µm	M=Internal 116/Flat with bayonet MN= Internal 116/Flat no bayonet CK= 126/Flat with 2 bayonet PK= Special Flange CM: Only for F25	S=Silicone E=EPDM N=Nitrile F=FKM

PP Melt Blown Filter Cartridge

MicroPure Classic Series

PP melt blown filter cartridge is made of 100% pure PP particles without any adhesive, and is made by fiber hot-melt entanglement. The unique three-layer gradient pore size design forms a three-dimensional slag filtering effect, which has the characteristics of high porosity, high retention efficiency, high dirt holding capacity and low pressure difference. The pure PP structure makes it have extensive chemical compatibility and is suitable for the filtration of strong acid, strong alkali and organic solvents. The advanced fiber hot-melt technology can meet the various needs of customers for this product.

Technical Parameters


	Filter Membrane	PP
	Cage/End Cap	PP
Materials of	Core	PP
Construction	Connection	DOE/222Flat/222Fin/226Fin/226Flat
	O-rings	S/F/N/E
	Surface Treatment	Flat/Orange/Embossing/Grooved
	Outer Diameter	63 / 114 mm
Filter Dimensions	Inner Diameter	28 / 30 mm
Difficultions	Removal Rating	0.5µm/1µm/5µm/10µm/20µm/30µm 50µm/75µm/100µm/150µm
	Length	10"/20"/ 30"/ 40"/ 50"/ 60"/ 70"
Operating Conditions	Operating Temperature and pH	<82°C, pH 1-13
Conditions	Maximum Differential Pressure (forward)	3.5 bar @ 20°C,1.0 bar @ 82°C
	Recommended Replacement Differential Pressure	2.4 bar @ 20°C

Features

- Gradient pore size structure, high dirt holding capacity
- $\bullet \ \ \mbox{High-purity PP material, extensive chemical resistance}$
- Uniform void distribution, high filtering effect, large area
- Double the service life and low replacement costs
- Inner ironing treatment and PP core structure to avoid fiber precipitation
- NFS certification



Typical Liquid Flow Rates



Typical Applications


- Desalination, RO reverse osmosis prefiltration
- Electronics industry: prefiltration of pure water and ultra-pure water
- Power plant water treatment: filtration of cooling water, boiler brine removal
- Pretreatment of biopharmaceutical, water for injection and other pure water
- Filtration of drinking water: filtration of food, beverages and bottled water
- Filtration of electroplating solutions
- Pretreatment of oil/gas industry, filtration of oilfield reinjection water

Ordering Information

	Length	Pore Size	Adaptor	Seal Material	Appearance	Core
CRMPP	010=10" 020=20" 030=30" 040=40" 050=50" 060=60" 070=70"	0050=0.5μm 0100=1μm 0500=5μm 1000=10μm 2000=20μm 3000=30μm 5000=50μm 10000=100μm	Blank=DOE 2=222/Flat 3=222/Fin 6=226/Fin 7=226/Flat	Blank=Non S=Silicone E=EPDM N=Nitrile F=FKM	B=Big and Fat (Conventional Surface) HG=Fluffy Surface (High Dirt Holding) G=Grooved Surface E=Embossing Surface O=Orange Surface Blank=Conventional Surface	P=PP Blank=Non

String wound filt liquid with low vi made of PP yarn skeleton accordin

Typical Liquid Flow Rates



Features

- Wide selection of filter materials and core materials ensure superior performance
- $\bullet \;$ Gradient pore structure effectively removes particles of different size
- Long service life and save operating costs
- Extensive chemical resistance
- High dirt holding capacity

String Wound Filter Cartridge

AquaPure Series

String wound filter cartridge is a kind of depth filter, which is used for the filtration of liquid with low viscosity and low impurity content. The string wound filter cartridge is made of PP yarn or absorbent cotton or glass fiber yarn precisely wound on the porous skeleton according to a specific density, which has a honeycomb structure with a sparse outside and a dense inside. It can effectively remove suspended matter, particles, rust and other impurities in the fluid, and has very excellent filtration performance.

Technical Parameters

	Filter Membrane	PP/Absorbent Cotton/Glass Fiber
	Cage/End Cap	PP
Materials of Construction	Core	PP/304/316 Stainless Steel
Construction	Connection	DOE/222Flat/222Fin/226Fin/226Flat
	O-rings	S/F/N/E
	Outer Diameter	63 / 114 mm
Filter	Inner Diameter	28 / 30 mm
Dimensions	Removal Rating	0.5µm/1µm/5µm/10µm/20µm/30µm 50µm/75µm/100µm/150µm
	Length	10"/20"/ 30"/ 40"/ 50"/ 60"/ 70"
	Operating Temperature	PP<82°C Absorbent Cotton/Glass Fiber<120°C
Operating Conditions	Maximum Differential Pressure (forward)	4.2 bar @ 20°C
	Recommended Replacement Differential Pressure	2.4 bar @ 20°C
	PH	1-13

Typical Applications

- Pretreatment of biopharmaceutical high-purity water/water for injection
- RO pretreatment
- Petrochemical/power plant filtration
- Filtration of industrial oils
- Filtration of microelectronics ultra-pure water
- Filtration of fine chemical acid and alkali or organic solutions

	Material	Length	Pore Size	Core	Adaptor	Seal Material
	PP	010=10"	005 = 0.5μm	P=PP Skeleton	Blank=DOE	Blank=Non
	CT=Absorbent Cotton	020=20"	010 = 1μm	S=304SS Skeleton	2=222/Flat	S=Silicone
CRW	GF=Glass Fiber	030=30" 040=40"	050 = 5μm 100 = 10μm	L=316L Skeleton	3=222/Fin	E=EPDM
		040=40	200 = 10μm		6=226/Fin	F=FKM
		060=60"	300 = 30μm		7=226/Flat	N=Nitrile
		070=70"	500 = 50μm			Rubber
			H10 = 100µm			

High Flow PP Pleated Filter Cartridge


600 FlowPure-HF Series

600 FlowPure-HF high flow series pleated filter cartridge has a large diameter of 6inch(152mm), no core, and a single opening design. The large diameter ensures a larger filtration area and increases the flux, thereby reducing the number of filters and cartridges. Long service life and high flow rate reduce investment costs.

Technical Parameters

	Filter Membrane	PP/Glass Fiber
	Support/Shell	PP/GF Reinforced PP
Materials of Construction	End Cap	Reinforced PP
	O-rings	S/E/N
	Welding Method	Hot-melt Welding(No Adhesive)
	Outer Diameter	152mm(6")
Filter Dimensions	Length	20"/40"/60"
	Connection	Carrying Handle Type
	Operating	Ordinary PP Shell <80℃
Operating	Temperature	GF Reinforced PP Shell <120℃
Conditions	Maximum Differential Pressure	3.4 bar @ 82°C
	Recommended Replacement Differential Pressure	2.4 bar @ 20°C

High Flow Dimensional Drawing





Features

- Filter material has gradient pore structure
- Flow rate is 5 times that of ordinary pleated filter cartridge
- 50% less investment in the overall system
- Uses hot-melt welding technology

Typical Applications

- RO prefiltration, desalination plants and industrial process water
- Filtration of power plants, condensed water
- Drinking water, food and beverage, edible oil, etc.
- Filtration of oilfield reinjection water
- High-purity water process pretreatment in biopharmaceutical industry

Ordering Information

	Filter Material	Length	Removal Rating	Seal Material	Appearance
CRHF	PP GF=Glass Fiber	020 = 20° 040 = 40° 060 = 60°	100 = 1μm 450 = 4.5μm 600 = 6μm 1000 = 10μm 2000 = 20μm 5000 = 50μm 7000 = 70μm 10000= 100μm	S=Silicone E=EPDM N=Nitrile Rubber	B=Bind Belt Type C=Shell Type



NMO Series

Nylon monofilament filter bag is made of filter material woven by a single nylon fiber, and the structure is welded and fused or wire stitched. It has good compressive strength, smooth surface and is not easy to deform, easy to clean and can be used repeatedly. It is suitable for the fields with low removal rating requirements such as coarse filtration or prefiltration as well as liquid filtration with high impurity content.

Technical Parameters

	Filter Membrane	Nylon Monofilament cloth	
Materials of Construction	Connection Material	PP/Polyester/Stainless Steel	
	Technique	Wire Stitching/Hot-melt	
	Bag 1	Φ7"×17" L; 0.25 m²	
Filter	Bag 2	Ф7″×32″ L; 0.50 m²	
Dimensions	Bag 3	Ф4″×9″ L; 0.09 m²	
	Bag 4	Ф4″×16″ L; 0.16 m²	
Operating Conditions	Operating Temperature	<160°C	

Features `

- Fixed aperture, constant value filtration
- High tensile strength
- Easy to clean and reusable
- Can intercept all kinds of hard particle impurities
- Suitable for high viscosity materials

Typical Applications

- Prefiltration of water treatment
- Petrochemical industry
- Light industry food and beverages
- Solid-liquid separation and treatment in metallurgy, electronics, pharmaceutical and other fields
- Water treatment for biopharmaceutical industry

	Material	Pore Size	Outer Filter Cloth	Dimension	Interface Ring Material
FB	NMO=Nylon Monofilament Cloth	025=25µm 050=50µm 075=75µm 100=100µm 200=200µm 300=300µm 400=400µm 600=600µm 800=800µm	P=Non N=Nylon Monofilament Cloth	01=Bag 1(Φ7"*17"L) 02=Bag 2(Φ7"*32"L) 03=Bag 3(Φ4"*9"L) 04=Bag 4(Φ4"*16"L)	P=PP Ring E=Polyester Ring S=Stainless Steel Ring

Needle Felt Filter Bag

PP/PE/PTFE Series

Needle felt filter bag is made of 100% ultra-fine pure fiber needle felt with three-dimensional deep filtration structure. It is highly fluffy, has a tortuous interior, and has super air permeability and high dirt holding capacity. It can effectively capture fine solid particles and colloidal particles, and is not easy to block. The surface of the felt cloth adopts high temperature heat treatment, singeing and calendering treatment to avoid pollution caused by fiber shedding, and to avoid the excessive blockage of the filter holes caused by the traditional rolling treatment and shorten the service life. The filter bag is available in different materials such as pure PP, polyester and PTFE to meet your different filtering

Technical Parameters

	Filter Membrane	PP/PE/PTFE Needle Felt
Materials of	Connection Material	PP/Polyester/Stainless Steel
Construction	Technique	Wire Stitching/Hot-melt
	Bag 1	Φ7″×17″ L; 0.25 ㎡
Filter	Bag 2	Ф7″×32″ L; 0.50 m²
Dimensions	Bag 3	Ф4″×9″ L; 0.09 m²
	Bag 4	Ф4"×16" L; 0.16 m²
Operating Conditions	Operating Temperature	PP < 90°C PE < 160°C PTFE < 240°C

Features

- High permeability and high dirt holding capacity
- No hair removal on the surface, no chemical release
- Hot-melt seamless welding
- Material does not contain silicone oil
- High flow rate, low pressure difference



Typical Applications

- Prefiltrration of water treatment
- · Food and beverage industry
- · Filtration of pretreatment precleaning, electrophoretic paint, spray water, circulating water, etc. in the automotive industry
- Solvent purification and treatment, resin and wastewater treatment, etc. in the paint industry
- Water treatment for biopharmaceutical industry

Ordering Information

	Material	Pore Size	Outer Filter Cloth	Dimension	Interface Ring Material
FB	PO=PP Needle Felt PE=Polyester Needle Felt PTFE	025 = 25µm 050 = 50µm 075 = 75µm 100 = 100µm 200 = 200µm 300 = 300µm 400 = 400µm 600 = 600µm 800 = 800µm	P=Non PM=PP Multifilament Cloth EM=Polyester Multifilament Cloth	01=Bag 1(Ф7"*17"L) 02=Bag 2(Ф7"*32"L) 03=Bag 3(Ф4"*9"L) 04=Bag 4(Ф4"*16"L)	P=PP Ring E=Polyester Ring S=Stainless Steel Ring

Depth Filter Sheets and Modules

CRDS Series



Typical Applications

- Clarification and filtration of raw material liquid
- Filtration of vaccines
- Decarbonization filtration of LVP
- Filtration of blood products
- Sterilization filtration of Chinese medicine liquid
- Fine chemical industry
- Food and beverages

CRDS series deep membrane stack filter cartridge compounded with high-purity lignocellulose and inorganic filter aids (diatomite, etc.). Its special three-dimensional space of channel makes each square meter filter area has thousands of square meters of internal surface area, and has a very high dirt holding capacity and filtration quality. Diatomite and resin can not only improve the permeability of the filter plate, but also improve the overall adsorption of the filter plate, and effectively filter the impurities and microorganisms in the liquid through mechanical interception, electrostatic adsorption, etc. The deep stacking technology effectively solves the problems of material and liquid loss, heavy workload and high investment cost of the traditional plate and frame filter.

Technical Parameters

	Filter Membrane	Cellulose/Diatomite/Resin,etc.
Materials of Construction	Skeleton/Support Separation Plate	PP
	O-rings	See ordering information selection table
Operating	Maximum Operating Temperature	80°C 1.0 bar/80°C
Conditions	Maximum Differential Pressure	2.0 bar/25°C
Product Guarantee	 The filter cartridges are a grade clean environment. Strictly implement ISO 90 management system All filter cartridge materials 	t 001:2015 quality

Features

- Thermal disinfection or hot filter fluid has no adverse effect on the filter board
- The operation is simple and reliable, and the solid outer skeleton design makes the filter cartridge not damaged during installation and disassembly
- High-purity cellulose greatly reduce the ion precipitation grade and endotoxin
- The three-dimensional space structure and the porous structure of diatomite greatly improve the dirt holding capacity and service life
- Flushable under certain conditions and extend the service life

	Pore Size	Connection	Outer Diameter	O-rings	Structural Unit
CRDS	002=0.2-0.4µm 004=0.4-0.6µm 006=0.6-1µm 010=1-3µm 020=2-5µm 030=3-7µm 100=10-15µm 250=25-30µm 400=40-50µm 500=50-60µm	D=Double Open End S=Single Open End (For 8"only)	08=8" 12=12" 16=16"	S=Silicone Rubber E=EPDM Rubber V=Fluorine Rubber	8=Layer 8 (8" DOE) 9=Layer 9 (12", 16") 12=Layer 12 (12", 16") 15=Layer 15 (12", 16") 16=Layer 16 (12", 16")

Metal pleated Filter Cartridge

CRPM Series

CRPM series metal pleated cartridge filter is pleated by special stainless steel wire mesh. The pleated process enables the filter cartridge to have a large filtration area, high dirt holding capacity and high flow rate. Each seal adopts the argon arc welding process without leakage. Under the high temperature and high pressure filtration environment, the stainless steel filter cartridge has excellent performance, and the filter cartridge can be repeatedly cleaned.

Technical Parameters

	Removal Rating	1 μm, 3 μm, 5 μm, 10 μm
	Length	10", 20", 30", 40"
Materials of	Filter Membrane	Stainless Steel Wire Mesh
Construction	Cage/End Cap/Core	Stainless Steel
	O-rings	See filter cartridge selection table
Operating	Maximum Operating Temperature	200°C
Conditions	Maximum Differential Pressure	5 bar



Typical Applications

- Filtration of steam
- · Filtration of oxidizing liquid
- Filtration of high viscosity liquid
- Filtration of liquid and gases under high temperature and high pressure condition
- Decarbonization filtration of material liquid
- Filtration of strong acid and alkali liquid

- · Pore distribution, good permeability
- No media shedding
- Strong corrosion resistance, not easy to damage
- Washable and long service life

Features

- Mental material has higher mechanical strength
- Wide temperature range that can be used

Ordering Information

	Material	Length	Removal Rating	End Cap	Seal Material
CRPM	4=304 Stainless Steel 6=316L Stainless Steel	005 = 5" 010 = 10" 020 = 20" 030 = 30" 040 = 40"	100 = 1μm 300 = 3μm 500 = 5μm H10 = 10μm	2=222/Flat 3=222/Fin 6=226/Fin 7=226/Flat	S=Silicone E=EPDM F=Fluorine N=Nitrile Rubber

Metal Titanium Rod Filter Cartridge

CRST Series

CRST series metal titanium rod filter cartridge is a microporous filter cartridge formed by high-purity titanium powder and sintered at high temperature and high vacuum. It has high porosity, good mechanical properties, high corrosion resistance, good chemical compatibility, no media shedding, very low leachables, can be repeatedly cleaned and reused, and low operating costs.

Technical Parameters

	Removal Rating	1μm/3μm/5μm/10μm
	Length	5"/10"/20"
Materials of Construction	Filter Membrane	Sintered Titanium Powder
construction.	Cage/End Cap/Core	High-purity Titanium
	O-rings	See filter cartridge selection table
Operating	Maximum Operating Temperature	≤280°C
Conditions	Maximum Differential Pressure	3.0 bar



Features

Эл. почта: info@tisys.ru info@tisys.kz info@tisys.by

- Tubular porous structure
- High mechanical strength
- Good temperature resistance
- Can be repeatedly cleaned and reused

Typical Applications

- Decarbonization filtration in concentration ration of LVP, small injection and oral liquid
- Impurity removal and filtration of API production, sterilization filtration of fermentation broth
- Filtration of corrosive reagents
- Filtration of high temperature liquid and steam

	Material	Length	Removal Rating	End Cap	Seal Material
CRS	T = Titanium	005 = 5" 010 = 10" 020 = 20"	100 = 1 μm 300 = 3 μm 500 = 5 μm H10 = 10 μm	2=222/Flat 7=226/Flat N=M30 Threaded Connection L=M20 Threaded Connection	S=Silicone E=EPDM F=Fluorine N=Nitrile Rubber